KR-800/RM-800

Auto Kerato-Refractometer Auto Refractometer





KR-800 / RM-800 Premium value based upon 60 years of experience

KR-800 & RM-800 Auto Kerato-Refractometers

The KR-800 & RM-800 Auto Kerato-Refractometers incorporate the very latest design technology and ergonomics. The units feature a bright 8.5 inch color touch screen panel to control the main functions and an improved joystick operation due to a 23% reduction in weight from the previous models.

Topcon systems have been renowned for their accuracy due to the proven Rotary Prism Technology. These units feature this technology as well, so accurate and stable measurements will be the norm. With over 60 years of experience, Topcon is always striving to deliver the very best in terms of instrumentation which reflects the best in design, functionality & accuracy.

Features

- Topcon's Rotary Prism Technology for unparalleled accuracy
- Simple to use
- Compact and modern design
- 8.5 inch LCD touch screen panel
- Connectable LAN




Simple to use
















Wide 8.5 inch touch screen control panel

The KR-800/RM-800 is equipped with a wide 8.5 inch color touch screen panel featuring clear, easy to read icons allowing the operator full control of the unit during the measuring process.

Smooth movements

Topcon's ability to engineer a weight reduction of approximately 23% has contributed to a smoother operation of the unit when using the X-Y control lever. This all contributes to a more fluid measuring process and a more accurate end result.

One touch lock

Both instruments are fitted with a quick and easy to use lock.

Easy load printer

The printer on the units saves you time by eliminating the need to load new printer rolls. You simply drop the new roll and print.

Unparalleled accuracy

The unique* Rotary Prism Technology, exclusively available at Topcon, allows for unparalleled precision and reliability. This system measures a wider area of the retina while permitting measurements through pupils as small as 2 mm in diameter.

*The off-centered ring target can cover areas that may be occluded by a small pupil. The angle of the prism enables a wider area of the retina to be measured.

Printed data

A complete print out of the acquired data is available at the touch of a button. The operator can choose among three printed data formats:

- **1. All:** Prints out all measurements (up to 10 readings per eye);
- **2. Ave:** Prints out an averaged value for all measurements;
- **3. Classic:** Outputs data of all (up to 10) measurements for both eyes and the results of one measurement and the average measured value of horizontal and vertical corneal curvature for both eyes.


The general addition and estimated near PD value can be printed out.

Specifications

| | | KR-800 | RM-800 |
|-------------------------------|-------------------------------------|--|--------|
| Refractive power measurement | Spherical refractive power | -25D to +22D (0.12D/0.25D steps) ** | |
| | Cylindrical refractive power | 0D to ±10D (0.12D/0.25D steps) ** | |
| | Astigmatic axial angle | 0° to 180° (1° / 5° steps) | |
| | Minimal measurable pupil diameter | φ 2 mm | |
| Corneal curvature measurement | Corneal curvature radius | 5.00 mm to 10.00 mm (0.01 mm step) | |
| | Corneal refractive power | 67.50D to 33.75D (0.12D/0.25D steps) (where, corneal refractive power = 1.3375) | |
| | Corneal astigmatic refractive power | 0D to ±10D (0.12D/0.25D steps) | |
| | Corneal astigmatic axial angle | 0° to 180° (1°/5° steps) | |
| PD measurement range | | 20 mm to 85 mm (0.5 mm step) | |
| Data transport terminal | | USB (import) / RS-232C (export) / LAN (export) | |
| Dimensions | | 317 - 341 mm (W) x 521 - 537 mm (D) x 447 - 477 mm (H) | |
| Weight | | 15 kg | |
| Power supply | | 100-240V AC, 50-60Hz, 30-70VA | |
| | | | |

※-25D ≤ spherical refractive power + cylindrical refractive power or spherical refractive power + cylindrical refractive power ≤ +22D

System chart



IMPORTANT

Subject to change in design and/or specifications without advanced notice. In order to obtain the best results with this instrument, please be sure to review all user instructions prior to operation.

Medical device Class I. Manufacturer: Topcon Corporation.









Topcon Furone Medical B.V.

Essebaan 11; 2908 LJ Capelle a/d LJssel; P.O. Box 145; 2900 AC Capelle a/d LJssel; The Netherlands Phone: +31-(0)10-4585077; Fax: +31-(0)10-4585045 E-mail: medical@topcon.eu; www.topcon-medical.eu

Topcon Danmark

Praestemarksvej 25; 4000 Roskilde, Denmark Phone: +45-46-327500; Fax: +45-46-327555 E-mail: info@topcon.dk www.topcon.dk

Topcon Scandinavia A.B.

Neongatan 2; P.O. Box 25; 43151 Mölndal, Sweden Phone: +46-(0)31-7109200; Fax: +46-(0)31-7109249 E-mail: medical@topcon.se; www.topcon.se

Topcon España S.A.

HEAD OFFICE; Frederic Mompou, 4; 08960 Sant Just Desvern; Barcelona, Spain Phone: +34-93-4734057; Fax: +34-93-4733932 E-mail: medica@topcon.es; www.topcon.es

Topcon Italy

Viale dell' Industria 60; 20037 Paderno Dugnano, (MI) Italy Phone: +39-02-9186671; Fax: +39-02-91081091 E-mail: info@topcon.it; www.topcon.it

Topcon France

BAT AI; 3 route de la révolte, 93206 Saint Denis Cedex Phone: +33-(0)1-49212323; Fax: +33-(0)1-49212324 E-mail: topcon@topcon.fr; www.topcon-medical.fr

Topcon Deutschland GmbH Hanns-Martin-Schleyer Strasse 41;

D-47877 Willich, Germany Phone: (+49) 2!54-885-0; Fax: (+49) 2!54-885-177 E-mail: info@topcon-medical.de; www.topcon-medical.de

Topcon Polska Sp. z o.o.

ul. Warszawska 23; 42-470 Siewierz; Poland Phone: +48-(0)32-670-50-45; Fax: +48-(0)32-671-34-05 www.topcon-polska.pl

Topcon Great Britain Ltd.

Topcon House; Kennet Side; Bone Lane; Newbury Berkshire RG14 5PX; United Kingdom Phone: +44-(0)1635-551120; Fax: +44-(0)1635-551170 E-mail: medical@topcon.co.uk, www.topcon.co.uk

Topcon Ireland

Unit 276, Blanchardstown; Corporate Park 2 Ballycoolin; Dublin 15, Ireland Phone: +353-18975900; Fax: +353-18293915 E-mail: medical@topcon.ie; www.topcon.ie